

Antibodies Product Datasheet

Product Name:	ID3 Antibody
Catalog Number:	C44184
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 68-117
Concentration:	1 mg/ml
Gene Name:	ID3 1R21 BHLHB25 HEIR1
Protein Name:	DNA-binding protein inhibitor ID-3 (Class B basic helix-loop-helix protein 25) (bHLHb25) (Helix-loop-helix protein HEIR-1) (ID-like protein inhibitor HLH 1R21) (Inhibitor of DNA binding 3)
Other Name:	
Human Swiss Prot No.:	Q02535
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ID3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year