

Antibodies Product Datasheet

Product Name:	IF2B1 Antibody
Catalog Number:	C44191
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	IGF2BP1 CRDBP VICKZ1 ZBP1
Protein Name:	Insulin-like growth factor 2 mRNA-binding protein 1 (IGF2 mRNA-binding protein 1) (IMP-1) (Coding region determinant-binding protein) (CRD-BP) (IGF-II mRNA-binding protein 1) (VICKZ family member 1) (Zip code-binding protein 1) (ZBP-1) (Zipcode-binding pr
Other Name:	
Human Swiss Prot No.:	Q9NZI8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IF2B1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year