

Antibodies Product Datasheet

Product Name:	IFI27 Antibody
Catalog Number:	C44195
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	IFI27
Protein Name:	Interferon alpha-inducible protein 27, mitochondrial (p27) (Interferon alpha-induced 11.5 kDa protein) (Interferon-stimulated gene 12a protein) (ISG12(a))
Other Name:	
Human Swiss Prot No.:	P40305
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IFI27 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year