

Antibodies Product Datasheet

Product Name:	IFI30 Antibody
Catalog Number:	C44196
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IFI30
Concentration:	1 mg/ml
Gene Name:	IFI30 GILT IP30
Protein Name:	IFI30
Other Name:	Gamma-interferon-inducible lysosomal thiol reductase (EC 1.8.-.-) (Gamma-interferon-inducible protein IP-30) (Legumaturain)
Human Swiss Prot No.:	P13284
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.06% sodium azide.
Specificity:	This antibody detects endogenous levels of IFI30 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year