

Antibodies Product Datasheet

Product Name:	IFIT2 Antibody
Catalog Number:	C44199
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 200-280
Concentration:	1 mg/ml
Gene Name:	IFIT2 CIG-42 G10P2 IFI54 ISG54
Protein Name:	Interferon-induced protein with tetratricopeptide repeats 2 (IFIT-2) (ISG-54 K) (Interferon-induced 54 kDa protein) (IFI-54K) (P54)
Other Name:	
Human Swiss Prot No.:	P09913
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IFIT2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year