

Antibodies Product Datasheet

Product Name:	IFN17 Antibody
Catalog Number:	C44206
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	IFNA17
Protein Name:	Interferon alpha-17 (IFN-alpha-17) (Interferon alpha-88) (Interferon alpha-I') (LeIF I) (Interferon alpha-T)
Other Name:	
Human Swiss Prot No.:	P01571
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IFN17 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year