

Antibodies Product Datasheet

Product Name:	IGDC4 Antibody
Catalog Number:	C44219
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IGDC4
Concentration:	1 mg/ml
Gene Name:	IGDCC4 DDM36 KIAA1628 NOPE
Protein Name:	Immunoglobulin superfamily DCC subclass member 4 (Neighbor of punc e11) (Protein DDM36) (hDDM36)
Other Name:	
Human Swiss Prot No.:	Q8TDY8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human IGDC4
Storage and Stability:	-20°C/1 year