

Antibodies Product Datasheet

Product Name:	IGFR1 Antibody
Catalog Number:	C44224
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IGFR1
Concentration:	1 mg/ml
Gene Name:	IGFLR1 TMEM149 U2AF1L4
Protein Name:	IGFR1
Other Name:	IGF-like family receptor 1 (Transmembrane protein 149) (U2 small nuclear RNA auxiliary factor 1-like 4)
Human Swiss Prot No.:	Q9H665
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.69% sodium azide.
Specificity:	This antibody detects endogenous levels of IGFR1 at Human/Mouse
Storage and Stability:	-20°C/1 year