

Antibodies Product Datasheet

Product Name:	IGHA1 Antibody
Catalog Number:	C44225
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IGHA1
Concentration:	1 mg/ml
Gene Name:	IGHA1
Protein Name:	IGHA1
Other Name:	Ig alpha-1 chain C region
Human Swiss Prot No.:	P01876
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.35% sodium azide.
Specificity:	This antibody detects endogenous levels of IGHA1 at Human
Storage and Stability:	-20°C/1 year