

Antibodies Product Datasheet

Product Name:	IGS11 Antibody
Catalog Number:	C44229
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	IGSF11 BTIGSF CXADRL1 VSIG3
Protein Name:	Immunoglobulin superfamily member 11 (IgSF11) (Brain and testis-specific immunoglobulin superfamily protein) (Bt-IGSF) (V-set and immunoglobulin domain-containing protein 3)
Other Name:	
Human Swiss Prot No.:	Q5DX21
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IGS11 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year