

Antibodies Product Datasheet

Product Name:	IGSF1 Antibody
Catalog Number:	C44230
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IGSF1
Concentration:	1 mg/ml
Gene Name:	IGSF1 IGDC1 KIAA0364 PGSF2
Protein Name:	Immunoglobulin superfamily member 1 (IgSF1) (Immunoglobulin-like domain-containing protein 1) (Inhibin-binding protein) (InhBP) (Pituitary gland-specific factor 2) (p120)
Other Name:	
Human Swiss Prot No.:	Q8N6C5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human IGSF1
Storage and Stability:	-20°C/1 year