

Antibodies Product Datasheet

Product Name:	IKBL1 Antibody
Catalog Number:	C44233
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 70-150
Concentration:	1 mg/ml
Gene Name:	NFKBIL1 IKBL
Protein Name:	NF-kappa-B inhibitor-like protein 1 (Inhibitor of kappa B-like protein) (I-kappa-B-like protein) (IkappaBL) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 1)
Other Name:	
Human Swiss Prot No.:	Q9UBC1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IKBL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year