

Antibodies Product Datasheet

Product Name:	IL1AP Antibody
Catalog Number:	C44235
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IL1AP
Concentration:	1 mg/ml
Gene Name:	IL1RAP C3orf13 IL1R3
Protein Name:	IL1AP
Other Name:	Interleukin-1 receptor accessory protein (IL-1 receptor accessory protein) (IL-1RAcP) (Interleukin-1 receptor 3) (IL-1R-3) (IL-1R3)
Human Swiss Prot No.:	Q9NPH3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.183% sodium azide.
Specificity:	This antibody detects endogenous levels of IL1AP at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year