

Antibodies Product Datasheet

Product Name:	IL36G Antibody
Catalog Number:	C44239
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	IL36G IL1E IL1F9 IL1H1 IL1RP2 UNQ2456/PRO5737
Protein Name:	Interleukin-36 gamma (IL-1-related protein 2) (IL-1RP2) (Interleukin-1 epsilon) (IL-1 epsilon) (Interleukin-1 family member 9) (IL-1F9) (Interleukin-1 homolog 1) (IL-1H1)
Other Name:	
Human Swiss Prot No.:	Q9NZH8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IL36G Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year