

Antibodies Product Datasheet

Product Name:	INGR2 Antibody
Catalog Number:	C44256
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 250-330
Concentration:	1 mg/ml
Gene Name:	IFNGR2 IFNGT1
Protein Name:	Interferon gamma receptor 2 (IFN-gamma receptor 2) (IFN-gamma-R2) (Interferon gamma receptor accessory factor 1) (AF-1) (Interferon gamma transducer 1)
Other Name:	
Human Swiss Prot No.:	P38484
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	INGR2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year