

Antibodies Product Datasheet

Product Name:	INP5E Antibody
Catalog Number:	C44257
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 490-570
Concentration:	1 mg/ml
Gene Name:	INPP5E
Protein Name:	72 kDa inositol polyphosphate 5-phosphatase (EC 3.1.3.36) (Phosphatidylinositol 4,5-bisphosphate 5-phosphatase) (Phosphatidylinositol polyphosphate 5-phosphatase type IV)
Other Name:	
Human Swiss Prot No.:	Q9NRR6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	INP5E Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year