

Antibodies Product Datasheet

Product Name:	INSL6 Antibody
Catalog Number:	C44260
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	INSL6 RIF1
Protein Name:	Insulin-like peptide INSL6 (Insulin-like peptide 6) (Relaxin/insulin-like factor 1) [Cleaved into: Insulin-like peptide INSL6 B chain; Insulin-like peptide INSL6 A chain]
Other Name:	
Human Swiss Prot No.:	Q9Y581
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	INSL6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year