

Antibodies Product Datasheet

Product Name:	INSM2 Antibody
Catalog Number:	C44262
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human INSM2
Concentration:	1 mg/ml
Gene Name:	INSM2 IA6 Nbla106
Protein Name:	INSM2
Other Name:	Insulinoma-associated protein 2 (Zinc finger protein IA-6)
Human Swiss Prot No.:	Q96T92
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.133% sodium azide.
Specificity:	This antibody detects endogenous levels of INSM2 at Human/Mouse
Storage and Stability:	-20°C/1 year