

Antibodies Product Datasheet

Product Name:	INSRR Antibody
Catalog Number:	C44263
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 80-160
Concentration:	1 mg/ml
Gene Name:	INSRR IRR
Protein Name:	Insulin receptor-related protein (IRR) (EC 2.7.10.1) (IR-related receptor) [Cleaved into: Insulin receptor-related protein alpha chain; Insulin receptor-related protein beta chain]
Other Name:	
Human Swiss Prot No.:	P14616
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	INSRR Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year