

Antibodies Product Datasheet

Product Name:	INT6 Antibody
Catalog Number:	C44266
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 550-630
Concentration:	1 mg/ml
Gene Name:	INTS6 DBI1 DDX26 DDX26A
Protein Name:	Integrator complex subunit 6 (Int6) (DBI-1) (Protein DDX26) (Protein deleted in cancer 1) (DICE1)
Other Name:	
Human Swiss Prot No.:	Q9UL03
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	INT6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year