

Antibodies Product Datasheet

Product Name:	IP6K1 Antibody
Catalog Number:	C44270
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	IP6K1 IHPK1 KIAA0263
Protein Name:	Inositol hexakisphosphate kinase 1 (InsP6 kinase 1) (EC 2.7.4.21) (Inositol hexaphosphate kinase 1)
Other Name:	
Human Swiss Prot No.:	Q92551
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IP6K1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year