

Antibodies Product Datasheet

Product Name:	IPO13 Antibody
Catalog Number:	C44271
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IPO13
Concentration:	1 mg/ml
Gene Name:	IPO13 KIAA0724 RANBP13
Protein Name:	IPO13
Other Name:	Importin-13 (Imp13) (Karyopherin-13) (Kap13) (Ran-binding protein 13) (RanBP13)
Human Swiss Prot No.:	O94829
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.263% sodium azide.
Specificity:	This antibody detects endogenous levels of IPO13 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year