

Antibodies Product Datasheet

Product Name:	IPO5 Antibody
Catalog Number:	C44272
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	IPO5 KPNB3 RANBP5
Protein Name:	Importin-5 (Imp5) (Importin subunit beta-3) (Karyopherin beta-3) (Ran-binding protein 5) (RanBP5)
Other Name:	
Human Swiss Prot No.:	O00410
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IPO5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year