

Antibodies Product Datasheet

Product Name:	IQEC1 Antibody
Catalog Number:	C44281
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IQEC1
Concentration:	1 mg/ml
Gene Name:	IQSEC1 ARFGEP100 BRAG2 KIAA0763
Protein Name:	IQEC1
Other Name:	IQ motif and SEC7 domain-containing protein 1 (ADP-ribosylation factors guanine nucleotide-exchange protein 100) (ADP-ribosylation factors guanine nucleotide-exchange protein 2) (Brefeldin-resistant Arf-GEF 2 protein)
Human Swiss Prot No.:	Q6DN90
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.264% sodium azide.
Specificity:	This antibody detects endogenous levels of IQEC1 at Human/Mouse
Storage and Stability:	-20°C/1 year