

Antibodies Product Datasheet

Product Name:	IREB2 Antibody
Catalog Number:	C44283
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IREB2
Concentration:	1 mg/ml
Gene Name:	IREB2
Protein Name:	IREB2
Other Name:	Iron-responsive element-binding protein 2 (IRE-BP 2) (Iron regulatory protein 2) (IRP2)
Human Swiss Prot No.:	P48200
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.265% sodium azide.
Specificity:	This antibody detects endogenous levels of IREB2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year