

Antibodies Product Datasheet

Product Name:	IRGM Antibody
Catalog Number:	C44284
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 101-150
Concentration:	1 mg/ml
Gene Name:	IRGM IFI1 IRGM1 LRG47
Protein Name:	Immunity-related GTPase family M protein (EC 3.6.5.-) (Immunity-related GTPase family M protein 1) (Interferon-inducible protein 1) (LPS-stimulated RAW 264.7 macrophage protein 47 homolog) (LRG-47)
Other Name:	
Human Swiss Prot No.:	A1A4Y4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IRGM Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year