

Antibodies Product Datasheet

Product Name:	IRPL1 Antibody
Catalog Number:	C44291
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 340-420
Concentration:	1 mg/ml
Gene Name:	IL1RAPL1 OPHN4
Protein Name:	Interleukin-1 receptor accessory protein-like 1 (IL-1-RAPL-1) (IL-1RAPL-1) (IL1RAPL-1) (Oligophrenin-4) (Three immunoglobulin domain-containing IL-1 receptor-related 2) (TIGIRR-2) (X-linked interleukin-1 receptor accessory protein-like 1)
Other Name:	
Human Swiss Prot No.:	Q9NZN1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IRPL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year