

Antibodies Product Datasheet

Product Name:	ISG15 Antibody
Catalog Number:	C44292
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ISG15 G1P2 UCRP
Protein Name:	Ubiquitin-like protein ISG15 (Interferon-induced 15 kDa protein) (Interferon-induced 17 kDa protein) (IP17) (Ubiquitin cross-reactive protein) (hUCRP)
Other Name:	
Human Swiss Prot No.:	P05161
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ISG15 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year