

Antibodies Product Datasheet

Product Name:	ITLN1 Antibody
Catalog Number:	C44305
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	ITLN1 INTL ITLN LFR UNQ640/PRO1270
Protein Name:	Intelectin-1 (ITLN-1) (Endothelial lectin HL-1) (Galactofuranose-binding lectin) (Intestinal lactoferrin receptor) (Omentin)
Other Name:	
Human Swiss Prot No.:	Q8WWA0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ITLN1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year