

Antibodies Product Datasheet

Product Name:	JIP4 Antibody
Catalog Number:	C44312
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 590-670
Concentration:	1 mg/ml
Gene Name:	SPAG9 HSS KIAA0516 MAPK8IP4 SYD1 HLC6
Protein Name:	C-Jun-amino-terminal kinase-interacting protein 4 (JIP-4) (JNK-interacting protein 4) (Cancer/testis antigen 89) (CT89) (Human lung cancer oncogene 6 protein) (HLC-6) (JNK-associated leucine-zipper protein) (JLP) (Mitogen-activated protein kinase 8-intera
Other Name:	
Human Swiss Prot No.:	O60271
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	JIP4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year