

Antibodies Product Datasheet

Product Name:	K0753 Antibody
Catalog Number:	C44317
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K0753
Concentration:	1 mg/ml
Gene Name:	KIAA0753
Protein Name:	K0753
Other Name:	Uncharacterized protein KIAA0753
Human Swiss Prot No.:	Q2KHM9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.274% sodium azide.
Specificity:	This antibody detects endogenous levels of K0753 at Human/Mouse
Storage and Stability:	-20°C/1 year