

Antibodies Product Datasheet

Product Name:	K121L Antibody
Catalog Number:	C44322
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K121L
Concentration:	1 mg/ml
Gene Name:	KIAA1211L C2orf55
Protein Name:	K121L
Other Name:	Uncharacterized protein KIAA1211-like
Human Swiss Prot No.:	Q6NV74
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.259% sodium azide.
Specificity:	This antibody detects endogenous levels of K121L at Human
Storage and Stability:	-20°C/1 year