

Antibodies Product Datasheet

Product Name:	K132L Antibody
Catalog Number:	C44324
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K132L
Concentration:	1 mg/ml
Gene Name:	KIAA1324L EIG121L
Protein Name:	UPF0577 protein KIAA1324-like (Estrogen-induced gene 121-like protein) (hEIG121L)
Other Name:	
Human Swiss Prot No.:	A8MWY0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human K132L
Storage and Stability:	-20°C/1 year