

Antibodies Product Datasheet

Product Name:	K1407 Antibody
Catalog Number:	C44325
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K1407
Concentration:	1 mg/ml
Gene Name:	KIAA1407
Protein Name:	Coiled-coil domain-containing protein KIAA1407
Other Name:	
Human Swiss Prot No.:	Q8NCU4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human K1407
Storage and Stability:	-20°C/1 year