

Antibodies Product Datasheet

Product Name:	K1751 Antibody
Catalog Number:	C44330
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K1751
Concentration:	1 mg/ml
Gene Name:	KIAA1751
Protein Name:	K1751
Other Name:	Uncharacterized protein KIAA1751
Human Swiss Prot No.:	Q9C0B2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.297% sodium azide.
Specificity:	This antibody detects endogenous levels of K1751 at Human/Mouse
Storage and Stability:	-20°C/1 year