

Antibodies Product Datasheet

Product Name:	K1908 Antibody
Catalog Number:	C44331
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K1908
Concentration:	1 mg/ml
Gene Name:	KIAA1908
Protein Name:	Uncharacterized protein KIAA1908
Other Name:	
Human Swiss Prot No.:	Q96PY0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human K1908
Storage and Stability:	-20°C/1 year