

Antibodies Product Datasheet

Product Name:	K1C27 Antibody
Catalog Number:	C44334
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K1C27
Concentration:	1 mg/ml
Gene Name:	KRT27 KRT25C
Protein Name:	Keratin, type I cytoskeletal 27 (Cytokeratin-27) (CK-27) (Keratin-25C) (K25C) (Keratin-27) (K27) (Type I inner root sheath-specific keratin-K25irs3)
Other Name:	
Human Swiss Prot No.:	Q7Z3Y8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human K1C27
Storage and Stability:	-20°C/1 year