

Antibodies Product Datasheet

Product Name:	K1H2 Antibody
Catalog Number:	C44337
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 220-300
Concentration:	1 mg/ml
Gene Name:	KRT32 HHA2 HKA2 KRTHA2
Protein Name:	Keratin, type I cuticular Ha2 (Hair keratin, type I Ha2) (Keratin-32) (K32)
Other Name:	
Human Swiss Prot No.:	Q14532
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	K1H2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year