

Antibodies Product Datasheet

Product Name:	K22O Antibody
Catalog Number:	C44339
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 220-300
Concentration:	1 mg/ml
Gene Name:	KRT76 KRT2B KRT2P
Protein Name:	Keratin, type II cytoskeletal 2 oral (Cytokeratin-2P) (CK-2P) (K2P) (Keratin-76) (K76) (Type-II keratin Kb9)
Other Name:	
Human Swiss Prot No.:	Q01546
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	K22O Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year