

Antibodies Product Datasheet

Product Name:	K2C6B Antibody
Catalog Number:	C44342
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K2C6B
Concentration:	1 mg/ml
Gene Name:	KRT6B K6B KRTL1
Protein Name:	K2C6B
Other Name:	Keratin, type II cytoskeletal 6B (Cytokeratin-6B) (CK-6B) (Keratin-6B) (K6B) (Type-II keratin Kb10)
Human Swiss Prot No.:	P04259
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.119% sodium azide.
Specificity:	This antibody detects endogenous levels of K2C6B at Human/Mouse
Storage and Stability:	-20°C/1 year