

Antibodies Product Datasheet

Product Name:	K2C6C Antibody
Catalog Number:	C44343
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	KRT6C KRT6E
Protein Name:	Keratin, type II cytoskeletal 6C (Cytokeratin-6C) (CK-6C) (Cytokeratin-6E) (CK-6E) (Keratin K6h) (Keratin-6C) (K6C) (Type-II keratin Kb12)
Other Name:	
Human Swiss Prot No.:	P48668
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	K2C6C Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year