

Antibodies Product Datasheet

Product Name:	KAD8 Antibody
Catalog Number:	C44345
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	AK8 C9orf98
Protein Name:	Adenylate kinase 8 (EC 2.7.4.3)
Other Name:	
Human Swiss Prot No.:	Q96MA6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KAD8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year