

Antibodies Product Datasheet

Product Name:	KAT6A Antibody
Catalog Number:	C44348
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	KAT6A MOZ MYST3 RUNXBP2 ZNF220
Protein Name:	Histone acetyltransferase KAT6A (EC 2.3.1.48) (MOZ, YBF2/SAS3, SAS2 and TIP60 protein 3) (MYST-3) (Monocytic leukemia zinc finger protein) (Runt-related transcription factor-binding protein 2) (Zinc finger protein 220)
Other Name:	
Human Swiss Prot No.:	Q92794
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KAT6A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year