

Antibodies Product Datasheet

Product Name:	KAZD1 Antibody
Catalog Number:	C44350
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KAZD1
Concentration:	1 mg/ml
Gene Name:	KAZALD1 FKSG28 FKSG40 UNQ2945/PRO21184
Protein Name:	Kazal-type serine protease inhibitor domain-containing protein 1
Other Name:	
Human Swiss Prot No.:	Q96I82
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KAZD1
Storage and Stability:	-20°C/1 year