

Antibodies Product Datasheet

Product Name:	KBRS1 Antibody
Catalog Number:	C44352
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 50-130
Concentration:	1 mg/ml
Gene Name:	NKIRAS1 KBRAS1
Protein Name:	NF-kappa-B inhibitor-interacting Ras-like protein 1 (I-kappa-B-interacting Ras-like protein 1) (Kappa B-Ras protein 1) (KappaB-Ras1)
Other Name:	
Human Swiss Prot No.:	Q9NYS0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KBRS1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year