

Antibodies Product Datasheet

Product Name:	KCNH6 Antibody
Catalog Number:	C44370
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 230-310
Concentration:	1 mg/ml
Gene Name:	KCNH6 ERG2
Protein Name:	Potassium voltage-gated channel subfamily H member 6 (Ether-a-go-go-related gene potassium channel 2) (ERG-2) (Eag-related protein 2) (Ether-a-go-go-related protein 2) (hERG-2) (hERG2) (Voltage-gated potassium channel subunit Kv11.2)
Other Name:	
Human Swiss Prot No.:	Q9H252
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KCNH6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year