

Antibodies Product Datasheet

Product Name:	KCNKG Antibody
Catalog Number:	C44373
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 230-310
Concentration:	1 mg/ml
Gene Name:	KCNK16 TALK1
Protein Name:	Potassium channel subfamily K member 16 (2P domain potassium channel Talk-1) (TWIK-related alkaline pH-activated K(+) channel 1) (TALK-1)
Other Name:	
Human Swiss Prot No.:	Q96T55
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KCNKG Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year