

Antibodies Product Datasheet

Product Name:	KCNS1 Antibody
Catalog Number:	C44376
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 470-550
Concentration:	1 mg/ml
Gene Name:	KCNS1
Protein Name:	Potassium voltage-gated channel subfamily S member 1 (Delayed-rectifier K(+) channel alpha subunit 1) (Voltage-gated potassium channel subunit Kv9.1)
Other Name:	
Human Swiss Prot No.:	Q96KK3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KCNS1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year