

Antibodies Product Datasheet

Product Name:	KDM3B Antibody
Catalog Number:	C44379
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 230-310
Concentration:	1 mg/ml
Gene Name:	KDM3B C5orf7 JHDM2B JMJD1B KIAA1082
Protein Name:	Lysine-specific demethylase 3B (EC 1.14.11.-) (JmjC domain-containing histone demethylation protein 2B) (Jumonji domain-containing protein 1B) (Nuclear protein 5qNCA)
Other Name:	
Human Swiss Prot No.:	Q7LBC6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KDM3B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year