

Antibodies Product Datasheet

Product Name:	KDM5A Antibody
Catalog Number:	C44381
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 720-800
Concentration:	1 mg/ml
Gene Name:	KDM5A JARID1A RBBP2 RBP2
Protein Name:	Lysine-specific demethylase 5A (EC 1.14.11.-) (Histone demethylase JARID1A) (Jumonji/ARID domain-containing protein 1A) (Retinoblastoma-binding protein 2) (RBBP-2)
Other Name:	
Human Swiss Prot No.:	P29375
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KDM5A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year